



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-Epidermal Growth Factor Receptor Mouse Monoclonal Antibody

Cat No.: KD-10343

Aliases:

EGFR; Epidermal Growth Factor Receptor; ERBB1; ERBB; Receptor Tyrosine-Protein Kinase ErbB-1; Erb-B2 Receptor Tyrosine Kinase 1; Proto-Oncogene C-ErbB-1; EC 2.7.10.1; HER1; Epidermal Growth Factor Receptor (Avian Erythroblastic Leukemia Viral (V-Erb-B) Oncogene Homolog); Erythroblastic Leukemia Viral (V-Erb-B) Oncogene Homolog (Avian); Avian Erythroblastic Leukemia Viral (V-Erb-B) Oncogene Homolog; Epidermal Growth Factor Receptor Tyrosine Kinase Domain; Cell Proliferation-Inducing Protein 61; Cell Growth Inhibiting Protein 40; EGFR VIII; EC 2.7.10; NISBD2; PIG61; MENA

Background:

UniProt Entry: [P00533](#);NCBI Gene Entry: [1956](#)

Application Information

Molecular Weight: Predicted, 134 kDa, observed, 175 kDa Conaley: Mouse monoclonal antibody Conne ID: 24GB5715 Species Reacnivy: Human Appctanions Tsed: Western Blotting (WB)

Immunogen

Recombinant protein of human EGFR

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB):1:2,500-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-epidermal growth factor receptor antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .