



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-eIF-2A protein kinase 2 Mouse Monoclonal Antibody

Cat No.: KD-10350

Aliases:

eIF-2A protein kinase 2; Protein kinase RNA-activated; PKR; p68 kinase; P1/eIF-2A protein kinase

Background:

UniProt Entry: [P19525](#);NCBI Gene Entry: [5610](#)

Application Information

Molecular Weight: Predicted, 62 kDa, observed, 68 kDa Conaley: Mouse monoclonal antibody Conne ID: 24GB5880 Spcies Reacnivy: Human Appctanions Tsed: Western Blotting (WB)

Immunogen

Recombinant protein of human EIF2AK2

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB):1:2,500-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-EIF2AK2 antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit . EIF2AK2, eukaryotic translation initiation factor 2 alpha kinase 2.



Western blotting analysis using anti-EIF2AK2 antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit .