



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-Fibroblast Growth Factor 2 Rabbit Monoclonal Antibody

Cat No.: KD-10379

Aliases:

FGF2; Fibroblast Growth Factor 2; FGFB; Fibroblast Growth Factor 2 (Basic); Heparin-Binding Growth Factor 2; HBGF-2; FGF-2; BFGF; Basic Fibroblast Growth Factor BFGF; Basic Fibroblast Growth Factor; Prostatropin

Background:

UniProt Entry: [P09038](#);NCBI Gene Entry: [2247](#)

Application Information

Molecular Weight: Predicted, 31 kDa, observed, 18-24kDa Clonality: Rabbit monoclonal antibody Clone ID: 24GB4420 Species Reactivity: Human Applications: Western Blotting (WB), Flow Cytometry (FCM)

Immunogen

Recombinant protein of human FGF2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-FGF2 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit . FGF2, fibroblast growth factor 2.



Western blotting analysis using anti-fibroblast growth factor 2 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:2,000. Green, isotype control; red, FGF2.