



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-GATA2/3 Rabbit Monoclonal Antibody

Cat No.: KD-10396

Aliases:

GATA2; GATA Binding Protein 2; NFE1B; Endothelial Transcription Factor GATA-2; GATA-Binding Protein 2; MONOMAC; IMD21; DCML; GATA3; GATA Binding Protein 3; HDR; Trans-Acting T-Cell-Specific Transcription Factor GATA-3; GATA-Binding Factor 3; GATA-Binding Protein 3; HDRS

Background:

UniProt Entry: [P23769/P23771](#);NCBI Gene Entry: [2624/2625](#)

Application Information

Molecular Weight: Predicted, 48,51 kDa; observed, 51 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB5730 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human GATA2/3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-GATA2/GATA3 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-GATA2/3 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .