



SunLong Biotech Co.,LTD  
Tel: 0086-571- 56623320 Fax:0086-571- 56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

## **[KD-Validated] Anti-Phospho-GATA3(S308) Rabbit Monoclonal Antibody**

Cat No.: KD-10397

### **Aliases:**

GATA3; GATA Binding Protein 3; HDR; Trans-Acting T-Cell-Specific Transcription Factor GATA-3; GATA-Binding Factor 3; GATA-Binding Protein 3; HDRS

### **Background:**

UniProt Entry: [P23771](#);NCBI Gene Entry: [2625](#)

### **Application Information**

Molecular Weight: Predicted, 48 kDa, observed, 48 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB2360 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

### **Immunogen**

A synthesized peptide derived from human Phospho-GATA3(S308)

### **Isotype**

Rabbit IgG

### **Storage Buffer**

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### **Storage**


Store at -20 °C for one year.

### Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

### Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-Phospho-GATA3(S308) antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-Phospho-GATA3(S308) antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .