



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-Growth differentiation factor 3 Rabbit Monoclonal Antibody

Cat No.: KD-10402

Aliases:

GDF3; Growth Differentiation Factor 3; Growth/Differentiation Factor 3; MCOPCB6;
MCOP7; GDF-3; KFS3

Background:

UniProt Entry: [Q9NR23](#);NCBI Gene Entry: [9573](#)

Application Information

Molecular Weight: Predicted, 41 kDa, observed, 16 kDa Clonality: Rabbit monoclonal
antibody Clone ID: 24GB6110 Species Reactivity: Human, Mouse, Rat Applications
Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human GDF3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-growth differentiation factor 3 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit .