



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-HNRNPK Rabbit Monoclonal Antibody

Cat No.: KD-10462

Aliases:

HNRNPK; Heterogeneous Nuclear Ribonucleoprotein K; TUNP; HNRPK; CSBP; Transformation Upregulated Nuclear Protein; Transformation Up-Regulated Nuclear Protein; DC-Stretch Binding Protein; HnRNP K; AUKS

Background:

UniProt Entry: [P61978](#);NCBI Gene Entry: [3190](#)

Application Information

Molecular Weight: Predicted, 51 kDa, observed, 60 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB4825 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB), Flow Cytometry (FCM)

Immunogen

A synthesized peptide derived from human hnRNP K

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Protocols


For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Flow cytometric analysis of HNRNPK expression in HepG2 cells using HNRNPK antibody 1:2,000. Green, isotype control; red, HNRNPK.



Western blotting analysis using anti-HNRNPK antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-HNRNPK antibody 1:10,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .