



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-Hsc70 Rabbit Monoclonal Antibody

Cat No.: KD-10477

Aliases:

HSPA8; Heat Shock Protein Family A (Hsp70) Member 8; HSC70; HSP73; HSPA10; HSC71; Lipopolysaccharide-Associated Protein 1; Heat Shock Cognate 71 KDa Protein; Heat Shock 70kDa Protein 8; LPS-Associated Protein 1; LAP-1; Epididymis Secretory Sperm Binding Protein Li 72p; N-Myristoyltransferase Inhibitor Protein 71; Constitutive Heat Shock Protein 70; Epididymis Luminal Protein 33; Heat Shock Cognate Protein 54; Heat Shock 70 KDa Protein 8; Heat Shock 70kd Protein 10; Heat Shock 70kD Protein 8; EC 3.6.4.10; HEL-S-72p; HEL-33; HSC54; HSP71; NIP71; LAP1

Background:

UniProt Entry: [P11142](#);NCBI Gene Entry: [3312](#)

Application Information

Molecular Weight: Predicted, 71 kDa, observed, 70 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB2675 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human Hsc70

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-Hsc70 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-Hsc70 antibody 1:5000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20000 respectively. Image was developed using FeQ™ ECL Substrate Kit .