



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-KSR1 Rabbit Monoclonal Antibody

Cat No.: KD-10538

Aliases:

KSR1; Kinase Suppressor Of Ras 1; RSU2; KSR; Kinase Suppressor Of Ras; EC 2.7.11.1

Background:

UniProt Entry: [Q8IVT5](#);NCBI Gene Entry: [8844](#)

Application Information

Molecular Weight: Predicted, 102 kDa, observed, 120 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB6290 Species Reactivity: Human, Rat Applications Tested: Western Blotting (WB), Flow Cytometry (FCM)

Immunogen

A synthesized peptide derived from human KSR1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Protocols


For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Flow cytometric analysis of KSR1 expression in HeLa cells using KSR1 antibody 1:2,000. Green, isotype control; red, KSR1.



Western blotting analysis using anti-KSR1 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-KSR1 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .