



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-Leukotriene A4 Hydrolase Mouse Monoclonal Antibody

Cat No.: KD-10563

Aliases:

Leukotriene A4 Hydrolase; Tripeptide Aminopeptidase LTA4H; Leukotriene A-4 Hydrolase; LTA-4 Hydrolase; Testicular Secretory Protein Li 27; Leukotriene A(4) Hydrolase; EC 3.4.11.4; EC 3.3.2.6; LTA4

Background:

UniProt Entry: [P09960](#);NCBI Gene Entry: [4048](#)

Application Information

Molecular Weight: Predicted, 69 kDa, observed, 69 kDa Conaley: Mouse monoclonal antibody Conne ID: 24GB5540 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

Immunogen

Recombinant protein of humann LTA4H

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB):1:2,500-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-leukotriene A4 hydrolase antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit .