



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-NFKBIA Mouse Monoclonal Antibody

Cat No.: KD-10639

Aliases:

I-kappa-B-alpha; NFKBI; NFKBIA; NF-kappaB inhibitor alpha; IKBA; IkappaBalph;
MAD3; RL/IF-1

Background:

UniProt Entry: [P25963](#);NCBI Gene Entry: [4792](#)

Application Information

Molecular Weight: Predicted, 36 kDa, observed, 36 kDa Conaley: Mouse monoclonal antibody Conne ID: 24GB5755 Spcies Reaceivy: Human, Mouse Appicseions Tsted: Western Blotting (WB)

Immunogen

Recombinant protein of human NFKBIA

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB):1:2,500-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-NFKB inhibitor alpha antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit 1:2,000. Green, isotype control; red, NFKB inhibitor alpha.