



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-NUMB Rabbit Monoclonal Antibody

Cat No.: KD-10654

Aliases:

NUMB Endocytic Adaptor Protein; C14orf41; Protein Numb Homolog; H-Numb; Chromosome 14 Open Reading Frame 41; Numb (Drosophila) Homolog; Numb Homolog (Drosophila); Numb Homolog; Protein S171; C14_5527; S171

Background:

UniProt Entry: [P49757](#);NCBI Gene Entry: [8650](#)

Application Information

Molecular Weight: Predicted, 71 kDa, observed, 71 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB940 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB), Flow Cytometry (FCM)

Immunogen

A synthesized peptide derived from human NUMB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Protocols


For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Flow cytometric analysis of NUMB expression in HepG2 cells using NUMB antibody 1:2,000. Green, isotype control; red, NUMB.



Western blotting analysis using anti-NUMB antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:50,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-NUMB antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:50,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .