



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-PVR Rabbit Monoclonal Antibody

Cat No.: KD-10753

Aliases:

PVR; PVR Cell Adhesion Molecule; Poliovirus Receptor; CD155; NECL5; HVED; PVS; Nectin-Like Protein 5; Necl-5; NECL-5; TAGE4; Nectin-Like 5; CD155 Antigen; Tage4

Background:

UniProt Entry: [P15151](#);NCBI Gene Entry: [5817](#)

Application Information

Molecular Weight: Predicted, 45 kDa, observed, 70 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB2575 Species Reactivity: Human Applctanions Tsed: Western Blotting (WB), Immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Poliovirus Receptor

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Protocols


For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Immunocytochemical staining of HT-1080 cells with PVR antibody 1:1,000. Nuclei were stained blue with DAPI; PVR was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Low. Scale bar: 20 µm.



Western blotting analysis using anti-PVR antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-PVR antibody 1:5000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20000 respectively. Image was developed using FeQ™ ECL Substrate Kit .