



SunLong Biotech Co.,LTD  
Tel: 0086-571- 56623320 Fax:0086-571- 56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

## **[KD-Validated] Anti-Retinoic Acid Receptor Alpha Rabbit Monoclonal Antibody**

Cat No.: KD-10777

### **Aliases:**

Retinoic Acid Receptor Alpha; RAR-Alpha; NR1B1; RARalpha; RAR; Nuclear Receptor Subfamily 1 Group B Member 1; Nucleophosmin-Retinoic Acid Receptor Alpha Fusion Protein NPM-RAR Long Form; Retinoic Acid Receptor, Alpha Polypeptide; Retinoic Acid Receptor, Alpha; PML-DDX5-RARA Fusion Protein; PML-DDX5-RARA Fusion

### **Background:**

UniProt Entry: [P10276](#);NCBI Gene Entry: [5914](#)

### **Application Information**

Molecular Weight: Predicted, 51 kDa, observed, 55 kDa Clonality: Rabbit monoclonal antibody Clone ID: 24GB1160 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

### **Immunogen**

A synthesized peptide derived from human Retinoic Acid Receptor alpha

### **Isotype**

Rabbit IgG

### **Storage Buffer**

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### Storage


Store at -20 °C for one year.

### Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

### Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-Retinoic acid receptor alpha antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-Retinoic Acid Receptor alpha antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQ™ ECL Substrate Kit .