



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-RAS p21 protein activator 1 Rabbit Monoclonal Antibody

Cat No.: KD-10779

Aliases:

RASA1; RAS P21 Protein Activator 1; P120GAP; GAP; Ras GTPase-Activating Protein 1; P120RASGAP; CM-AVM; P120; RASA; RAS P21 Protein Activator (GTPase Activating Protein) 1; P120 RAS GTPase Activating Protein; Capillary Malformation-Arteriovenous Malformation; Triphosphatase-Activating Protein; GTPase-Activating Protein; Ras P21 Protein Activator; CMAVM1; RASGAP; RasGAP; CMAVM; PKWS

Background:

UniProt Entry: [P20936](#);NCBI Gene Entry: [5921](#)

Application Information

Molecular Weight: Predicted, 116 kDa, observed, 116 kDa Clonality: Rabbit monoclonal antibody Clone ID: 23GB4125 Species Reactivity: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human RASA1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-RAS p21 protein activator 1 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-RAS p21 protein activator 1 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .