



SunLong Biotech Co.,LTD  
Tel: 0086-571- 56623320 Fax:0086-571- 56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

---

## **[KD-Validated] Anti-PCNA Mouse Monoclonal Antibody**

Cat No.: KD-10957

### **Aliases:**

PCNA; Proliferating Cell Nuclear Antigen 2; Cyclin; DNA Polymerase Delta Auxiliary Protein; ATLD2

### **Background:**

UniProt Entry: [P12004](#);NCBI Gene Entry: [5111](#)

### **Application Information**

Molecular Weight: Predicted, 29 kDa, observed, 35 kDa Conaley: Mouse monoclonal antibody Conne ID: 24GB2805 Species Reactivy: Human, Mouse, Rat Applications Tested: Western Blotting (WB)

### **Immunogen**

Recombinant protein of human PCNA

### **Isotype**

Mouse IgG1 kappa

### **Storage Buffer**

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### **Storage**


Store at -20 °C for one year.

### Recommended Dilutions


Western Blotting (WB): 1:1,000-1:5,000

### Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-Proliferating cell nuclear antigen antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-Proliferating cell nuclear antigen antibody 1:5,000 and HRP-conjugated goat anti-mouse secondary antibody 1:20,000 respectively. Image was developed using FeQ™ ECL Substrate Kit .