



Goat Anti-Rat IgG H&L / Biotin antibody

SL0293G-Bio

Product Name	Goat Anti-Rat IgG H&L / Biotin
Chinese Name	生物素标记的羊抗大鼠 IgG H&L
Alias	Goat Anti-Rat IgG H&L (Biotin); Immunoglobulin G;
Immunogen Species	Goat
Clonality	Polyclonal
React Species	Rat, WB=1:2000-20000,IHC-P=1:200-1000,IHC-F=1:200-1000,ELISA=1:2000-20000,ICC/IF=1:10
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native Rat IgG
Lsotype	IgG
Purification	affinity purified by Protein G
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or d Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels betw important for our defence against microorganisms and the molecules are produced by B lymphocy The IgG molecule has two separate functions; to bind to the pathogen that elicited the response a destroy the antigen. The variability of the IgG pool is generated by somatic recombination and th given time point is estimated to be 1011 variants.
Product Detail	