

## Donkey Anti-Rabbit IgG H&L / AF350 antibody

SL0295D-AF350

**Product Name** Donkey Anti-Rabbit IgG H&L / AF350  
**Chinese Name** Alexa Fluor 350 标记的驴抗兔 IgG H&L  
**Alias** Donkey Anti-Rabbit IgG H&L (Alexa Fluor® 350); Immunoglobulin G;

**Specific References (1)** | SL0295D-AF350 has been referenced in 1 publications.



**[IF=9.685]** Song, Chen. et al. The phosphorylation of PHF5A by TrkA-ERK1/2-ABL1 cascade regulates centrosome separation. CELL DEATH DIS. 2023 Feb;14(2):1-14 **IHC,IF ; Mouse,Human.**

PubMed:36759599

**Immunogen Species** Donkey  
**Clonality** Polyclonal  
**React Species** Rabbit,  
 IF=1:100-1000,ICC/IF=1:100-1000,Flow-Cyt=1:100-1000  
**Applications** not yet tested in other applications.  
 optimal dilutions/concentrations should be determined by the end user.  
**Form** Liquid  
**Concentration** 2.0 mg/ml  
**immunogen** Native rabbit IgG  
**Lsotype** IgG  
**Purification** affinity purified by Protein G  
**Buffer Solution** 10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.  
**Storage** Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.  
**Product Detail** Immunoglobulin G (IgG), is one of the most abundant proteins in serum with



normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10<sup>11</sup> variants.