

Goat Anti-Rabbit IgG H&L / Streptavidin antibody

SL0295G-SA

Product Name Goat Anti-Rabbit IgG H&L / Streptavidin

Chinese Name 链霉亲和素标记山羊抗兔 IgG H&L

Alias Rabbit Immunoglobulin G; Rb IgG;

Specific References (1) | SL0295G-SA has been referenced in 1 publications.

[IF=3.391] Hui Yang. et al. A combined proteomic and metabolomic analyses of the priming phase during rat liver regeneration. Arch Biochem Biophys. 2020 Oct;693:108567 **WB ; Rat.**



PubMed:32898568

Immunogen Species Goat

Clonality Polyclonal

React Species Rabbit,

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Form Liquid

Concentration 1.0 mg/ml

immunogen Native Rabbit IgG

Lsotype IgG

Purification Affinity purified by Protein G

Buffer Solution 10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.

Storage Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Product Detail