



Mouse Anti-Rabbit IgG H&L / AF594 antibody

SL0295M-AF594

Product Name	Mouse Anti-Rabbit IgG H&L / AF594
Chinese Name	Alexa Fluor 594 标记的小鼠抗兔 IgG H&L
Alias	Mouse Anti-Rabbit IgG H&L (Alexa Fluor® 594); Immunoglobulin G;
Immunogen Species	Mouse
Clonality	Polyclonal
React Species	Rabbit, IF=1:100-1000,ICC/IF=1:100-1000,Flow-Cyt=1:100-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native rabbit IgG
Lsotype	IgG
Purification	affinity purified by Protein G
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	