

Mouse Anti-Rabbit IgG H&L / Gold antibody

SL0295M-Gold

Product Name	Mouse Anti-Rabbit IgG H&L / Gold
Chinese Name	胶体金标记的小鼠抗兔 IgG H&L
Alias	Mouse Anti-Rabbit IgG H&L (Gold); Immunoglobulin G;
Immunogen Species	Mouse
Clonality	Polyclonal
React Species	Rabbit, ICA=1:20-200,IEM=1:20-200
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	0.4 mg/ml
immunogen	Native rabbit IgG
Lsotype	IgG
Purification	affinity purified by Protein G
Buffer Solution	20 mM TBS (pH=7.4) with 1% BSA and 3% Proclin300.
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Product Detail	