



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

Goat Anti-Human IgG H&L whole serum antibody

SL0297Gs

Product Name	Goat Anti-Human IgG H&L whole serum
Chinese Name	羊抗人 IgG H&L 抗血清
Alias	Human Immunoglobulin G; hu IgG;
Immunogen Species	Goat
Clonality	Polyclonal
React Species	Human, WB=1:10000-50000 ELISA=1:10000-100000 IHC-P=1:10000-50000 IHC-F=1:10000-50000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Lyophilized powder
immunogen	Native Human IgG
Lsotype	IgG
Purification	Unpurified
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Product Detail	