



## Rabbit Anti-Human IgG H&L / AF594 antibody

SL0297R-AF594

<b>Product Name</b>	Rabbit Anti-Human IgG H&L / AF594
<b>Chinese Name</b>	Alexa Fluor 594 标记的兔抗人 IgG H&L
<b>Alias</b>	Rabbit Anti-Human IgG H&L (Alexa Fluor® 594); Immunoglobulin G;
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	Human, IF=1:100-1000
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Form</b>	Liquid
<b>Concentration</b>	2.0 mg/ml
<b>immunogen</b>	Native Human IgG
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
<b>Product Detail</b>	