

## Rabbit Anti-Pig IgG H&L / FITC antibody

SL0309R-FITC

**Product Name** Rabbit Anti-Pig IgG H&L / FITC  
**Chinese Name** FITC 标记的兔抗猪 IgG H&L  
**Alias** Rabbit Anti-Pig IgG H&L (FITC); Immunoglobulin G;

**Specific References (2)** | SL0309R-FITC has been referenced in 2 publications.

**[IF=2.791]** Zhang H et al. Construction and immunological evaluation of recombinant Newcastle disease virus vaccines expressing highly pathogenic porcine reproductive and respiratory syndrome virus GP3/GP5 proteins in pigs. *Veterinary Microbiology*,2019 108490. **ICF ; Pig.**



PubMed:doi:10.1016/j.vetmic.2019.108490

**[IF=2.524]** Liu et al. Resolution of the cellular proteome of the nucleocapsid protein from a highly pathogenic isolate of porcine reproductive and respiratory syndrome virus identifies PARP-1 as a cellular target whose interaction is critical for virus biology. (2015) *Vet.Microbiol.* 176:109-19 **IF ; Pig.**

PubMed:25614100

**Immunogen Species** Rabbit  
**Clonality** Polyclonal  
**React Species** Pig,  
 IF=1:100-1000,Flow-Cyt=1:100-1000,ICC/IF=1:100-1000  
**Applications** not yet tested in other applications.  
 optimal dilutions/concentrations should be determined by the end user.  
**Form** Liquid  
**Concentration** 2.0 mg/ml  
**immunogen** Native Pig IgG



---

<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	<p>This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.</p> <p>Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10<sup>11</sup> variants.</p>
<b>Product Detail</b>	