

Goat Anti-Chicken IgG H&L / Biotin antibody

SL0310G-Bio

Product Name Goat Anti-Chicken IgG H&L / Biotin

Chinese Name 生物素标记的羊抗鸡 IgG H&L

Alias Goat Anti-Chicken IgG H&L (Biotin); Immunoglobulin G;

Specific References (1) | SL0310G-Bio has been referenced in 1 publications.



[IF=4.171] Jiajun Yang. et al. Selenium enriched Bacillus subtilis yb-1114246 activated the TL chicken intestinal β -defensin 1 expression. 2021 May 18 **IHC ; Chicken.**

PubMed:34028482

Immunogen Species

Goat

Clonality

Polyclonal

React Species

Chicken,

WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000

Applications

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

Form

Liquid

Concentration

2.0 mg/ml

immunogen

Native Chicken IgG

Lsotype

IgG

Purification

affinity purified by Protein G

Buffer Solution

10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Attention

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic. Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 1-2 g/dl. It is important for our defence against microorganisms and the molecules are produced by B lymphocytes in response to an antigen. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to destroy the antigen. The variability of the IgG pool is generated by somatic recombination. The number of individual at a given time point is estimated to be 1011 variants.

Product Detail



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com
