



## Rabbit Anti-Chicken IgG H&L whole serum antibody

SL0310Rs

|                          |   |
|--------------------------|---|
| <b>Product Name</b>      | Rabbit Anti-Chicken IgG H&L whole serum   |
| <b>Chinese Name</b>      | 兔抗鸡 IgG H&L 抗血清   |
| <b>Alias</b>             | Chicken Immunoglobulin G  |
| <b>Immunogen Species</b> | Rabbit  |
| <b>Clonality</b>         | Polyclonal  |
| <b>React Species</b>     | Chicken,<br>WB=1:10000-50000 ELISA=1:10000-100000 IHC-P=1:10000-50000 IHC-F=1:10000-50000                                 |
| <b>Applications</b>      | not yet tested in other applications.<br>optimal dilutions/concentrations should be determined by the end user.           |
| <b>Form</b>              | Lyophilized powder  |
| <b>immunogen</b>         | Native Chicken IgG  |
| <b>Lsotype</b>           | IgG   |
| <b>Purification</b>      | Unpurified  |
| <b>Storage</b>           | Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.                                |
| <b>Attention</b>         | This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. |
| <b>Product Detail</b>    |   |