

Rabbit Anti-Avidin / Cy5 antibody

SL0312R-Cy5

Product Name	Rabbit Anti-Avidin / Cy5
Chinese Name	Cy5 标记的兔抗亲和素
Alias	Rabbit Anti-Avidin IgG H&L (Cy5)
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Avidin
Applications	IF=1:100-1000,Flow-Cyt=1:100-1000,ICC/IF=1:100-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native Avidin
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).
Product Detail	