

## Rabbit Anti-Avidin / Gold antibody

SL0312R-Gold

|                          |   |
|--------------------------|---|
| <b>Product Name</b>      | Rabbit Anti-Avidin / Gold   |
| <b>Chinese Name</b>      | 胶体金标记的兔抗亲和素   |
| <b>Alias</b>             | Rabbit Anti-Avidin IgG H&L (Gold)   |
| <b>Immunogen Species</b> | Rabbit  |
| <b>Clonality</b>         | Polyclonal  |
| <b>React Species</b>     | Avidin  |
| <b>Applications</b>      | ICA=1:20-200,IEM=1:20-200<br>not yet tested in other applications.<br>optimal dilutions/concentrations should be determined by the end user.  |
| <b>Form</b>              | Liquid  |
| <b>Concentration</b>     | 0.4 mg/ml   |
| <b>immunogen</b>         | Native Avidin   |
| <b>Lsotype</b>           | IgG   |
| <b>Purification</b>      | affinity purified by Protein A  |
| <b>Buffer Solution</b>   | 10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.   |
| <b>Storage</b>           | Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.   |
| <b>Attention</b>         | This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.<br>The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin). |
| <b>Product Detail</b>    |   |