



Rabbit Anti-Avidin / HRP antibody

SL0312R-HRP

Product Name Rabbit Anti-Avidin / HRP

Chinese Name 辣根过氧化物酶标记的兔抗亲和素

Alias Rabbit Anti-Avidin IgG H&L (HRP)

Immunogen Species Rabbit

Clonality Polyclonal

WB=1:1000-10000 ELISA=1:1000-10000 IHC-P=1:100-1000 IHC-F=1:100-1000

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Form Liquid

Concentration 2.0 mg/ml

immunogen Native Avidin

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.

Storage Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Product Detail The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).