

## Rabbit Anti-Avidin / PE-Cy7 antibody

SL0312R-PE-Cy7

<b>Product Name</b>	Rabbit Anti-Avidin / PE-Cy7
<b>Chinese Name</b>	PE-Cy7 标记的兔抗亲和素
<b>Alias</b>	Rabbit Anti-Avidin IgG H&L (PE-Cy7)
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Applications</b>	Flow-Cyt=1:100-1000 IF=1:100-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Form</b>	Liquid
<b>Concentration</b>	2.0 mg/ml
<b>immunogen</b>	Native Avidin
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).
<b>Product Detail</b>	