



Rabbit Anti-Avidin / RBITC antibody

SL0312R-RBITC

Product Name	Rabbit Anti-Avidin / RBITC
Chinese Name	罗丹明标记的兔抗亲和素
Alias	Rabbit Anti-Avidin IgG H&L (RBITC)
Immunogen Species	Rabbit
Clonality	Polyclonal IF=1:100-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native Avidin
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).