

Rabbit Anti-E1 glycoprotein antibody

SL0319R

Product Name	E1 glycoprotein
Chinese Name	风疹病毒 glycoprotein 抗体
Alias	E1; E1 envelope glycoprotein; Glycoprotein E1; Spike glycoprotein E1; POLS_RUBVM.
Research Area	immunology Bacteria and viruses
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Rubella virus) IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000 (Paraffin sections need antigen repair)
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	52kDa
Form	Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human Glycoprotein: 801-900/1063
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:Rubella virus)1M TBS(pH7.4) with 1% BSA, (predicted:Rubella virus)3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	