



Rabbit Anti-Bovine IgM whole serum antibody

SL0327Rs

Product Name	Rabbit Anti-Bovine IgM whole serum
Chinese Name	兔抗牛 IgM 抗血清
Alias	Bovine Immunoglobulin G
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Cow, WB=1:10000-50000 ELISA=1:10000-100000 IHC-P=1:10000-50000 IHC-F=1:10000-50000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Lyophilized powder
immunogen	Native Bovine IgM
Lsotype	IgG
Purification	Unpurified
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. IgM normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and is largely confined to plasma due to it's large size. Monomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma cells. Due to it's high valency IgM is more efficient than other isotypes is binding antigens with repeating epitopes (virus particles and red blood cells) and is more efficient than IgG in activating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.
Product Detail	