

Rabbit Anti-Hepatitis B Virus Surface Antigen antibody

SL0331R

Product Name	Hepatitis B Virus Surface Antigen
Chinese Name	肝炎病毒 S 蛋白抗体
Alias	HBSAG_HBVE1; Large envelope protein; L glycoprotein; L-HBsAg; LHB; Large S protein; Large surface protein; Major surface antigen.
Research Area	Bacteria and viruses
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Hepatitis B virus (HBVE1)) IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500 (Paraffin sections need antigen repair)
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	44kDa
Cellular localization	cytoplasmic The cell membrane
Form	Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from HBVE1 Hepatitis B Virus Surface Antigen
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:Hepatitis B virus (HBVE1))1M TBS(pH7.4) with 1% BSA, (predicted:Hepatitis B virus (HBVE1))3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Product Detail