

Rabbit Anti-digoxin whole serum antibody

SL0356Rs

Product Name	Rabbit Anti-digoxin whole serum
Chinese Name	地高辛抗血清
Alias	Digoxigenin; Lanoxin. 兔抗地高辛抗血清
Product Type	Small molecule anti
Immunogen Species	Rabbit
Clonality	Polyclonal
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	0.765kDa
Form	Liquid
immunogen	BSA conjugated digoxin
Lsotype	whole antiserum
Purification	Unpurified
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).
Product Detail	