

Goat Anti-Guinea Pig IgG H&L / Biotin antibody

SL0358G-Bio

Product Name Goat Anti-Guinea Pig IgG H&L / Biotin

Chinese Name 生物素标记的羊抗豚鼠 IgG H&L

Alias Goat Anti-Guinea Pig IgG H&L (Biotin); Immunoglobulin G;

Specific References (1) | SL0358G-Bio has been referenced in 1 publications.



[IF=4.362] JunXia Gao. et al. GLT-1 Knockdown Inhibits Ceftriaxone-Mediated Improvements in APP/PS1 AD Mice. Front Aging Neurosci. 2020; 12: 580772 **IHC ; M**

PubMed:33132901

Immunogen Species

Goat

Clonality

Polyclonal

React Species

Guinea Pig,

WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000

Applications

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

Form

Liquid

Concentration

2.0 mg/ml

immunogen

Native Guinea Pig IgG

Lsotype

IgG

Purification

affinity purified by Protein G

Buffer Solution

10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Attention

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic.

Product Detail

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 1-2 g/L. It is important for our defence against microorganisms and the molecules are produced by B lymphocytes in response to an antigen. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to destroy the antigen. The variability of the IgG pool is generated by somatic recombination of the variable regions of the heavy and light chain genes.



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

individual at a given time point is estimated to be 1011 variants.