



Rabbit Anti-FITC / Cy3 antibody

SL0366R-Cy3

Product Name	Rabbit Anti-FITC / Cy3
Chinese Name	Cy3 标记兔抗 FITC
Alias	fluorescein isothiocyanate; Fluorescein 5-isothiocyanate. 荧光素 Cy3 标记兔抗 FITC
Immunogen Species	Rabbit
Clonality	Polyclonal
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	1.0 mg/ml
Lsotype	
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	