

Mouse Anti-Rabbit IgM / PE-Cy7 antibody

SL0369M-PE-Cy7

Product Name	Mouse Anti-Rabbit IgM / PE-Cy7
Chinese Name	PE-Cy7 标记的小鼠抗兔 IgM
Alias	Mouse Anti-Rabbit IgM (PE-Cy7); Immunoglobulin M;
Immunogen Species	Mouse
Clonality	Polyclonal
React Species	Rabbit, Flow-Cyt=1:100-1000 IF=1:100-1000
Applications	IF=1:100-1000,Flow-Cyt=1:100-1000,ICC/IF=1:100-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native Rabbit IgM
Isotype	IgG
Purification	affinity purified by Protein G
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. Immunoglobulin M (IgM) normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and is largely confined to plasma due to its large size. Monomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma cells. Due to its high valency IgM is more efficient than other isotypes is binding antigens with repeating epitopes (virus particles and red blood cells) and is more efficient than IgG in activating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.
Product Detail	