

Rabbit Anti-Dog IgM / AF488 antibody

SL0373R-AF488

Product Name	Rabbit Anti-Dog IgM / AF488
Chinese Name	Alexa Fluor 488 标记的兔抗犬 IgM 抗体
Alias	Rabbit Anti-Dog IgM (Alexa Fluor® 488); Immunoglobulin M;
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Dog, IF=1:100-1000,Flow-Cyt=1:100-1000,ICC/IF=1:100-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	2.0 mg/ml
immunogen	Native Dog IgM
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	