



Rabbit Anti-Dog IgM / Gold antibody

SL0373R-Gold

Product Name	Rabbit Anti-Dog IgM / Gold
Chinese Name	胶体金标记的兔抗犬 IgM 抗体
Alias	Rabbit Anti-Dog IgM (Gold); Immunoglobulin M;
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Dog, ICA=1:20-200,IEM=1:20-200
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	0.4 mg/ml
immunogen	Native Dog IgM
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	