



Rabbit Anti-Dog IgM / HRP antibody

SL0373R-HRP

Product Name Rabbit Anti-Dog IgM / HRP

Chinese Name 辣根过氧化物酶标记的兔抗犬 IgM 抗体

Alias Rabbit Anti-Dog IgM (HRP); Immunoglobulin M;

Immunogen Species Rabbit

Clonality Polyclonal

React Species Dog,
WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,ELISA=1:1000-10000

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Form Liquid

Concentration 2.0 mg/ml

immunogen Native Dog IgM

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.

Storage Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Product Detail