

Rabbit Anti-BADH2 antibody

SL0612R

Product Name BADH2

Chinese Name BADH2(水稻)抗体

Alias Betaine-aldehyde dehydrogenase.

Product Type Botany 抗体

Research Area Cell biology

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted:Rice)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 55/172kDa

Cellular localization cytoplasmic

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from Rice BADH2: 101-200/503

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution (predicted:Rice)1M TBS(pH7.4) with 1% BSA, (predicted:Rice)3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Detail