



Rabbit Anti-Cobra poison Protein antibody

SL0617R

Product Name	Cobra poison Protein
Chinese Name	眼镜蛇蛇毒蛋白抗体
Research Area	immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Cobra) WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000 (Paraffin sections need antigen repair)
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	1mg/ml
immunogen	Cobra poison Protein(purified)
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:Cobra)1M TBS(pH7.4) with 1% BSA, (predicted:Cobra)3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	抗蛇毒抗体是由眼镜蛇毒液中采集、分离的原毒液血清免疫动物制备的抗体