

Rabbit Anti-Centromeric histone H3-like protein-1/Cy5 Conjugated antibody

SL0620R-Cy5

Product Name	Anti-Centromeric histone H3-like protein-1/Cy5
Chinese Name	Cy5 标记的组蛋白 H3 样蛋白抗体
Alias	Centromeric histone H3-like protein-1; centromere-specific H3 variant protein; A4PIF0_BRAJU;
Product Type	Botany 抗体
Research Area	Tumour immunology Developmental biology transcriptional regulatory factor Botany
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Brassica juncea) IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	13.6kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from Brassica juncea Centromeric histone H3-like protein-1 C-terminus
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail	Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com
