

Rabbit Anti-E.coli O6 antibody

SL0658R

Product Name	E.coli O6
Chinese Name	大肠杆菌埃希杆菌 O6 抗体
Alias	Escherichia coli O6; EPEC O6. 肠致病性大肠杆菌 O6
Research Area	immunology Bacteria and viruses
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Escherichia Coli) ELISA=1:5000-10000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	1mg/ml
immunogen	Escherichia coli O6 whole cell protein
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed This antibody eacts with Escherichia coli (E.coli), Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriae - bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, especially where there is an immature immune system or immunosuppression, or where the subtype of organism has acquired the ability to produce pathogenic toxins.
Product Detail	SWISS: N/A



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Gene ID:
N/A