



Rabbit Anti-Os05g0477200 protein antibody

SL0807R

Product Name Os05g0477200 protein

Chinese Name 水稻白叶枯病菌 Os05g0477200 蛋白抗体

Alias Oryza sativa subsp. japonica (Rice) Os05g0477200.

Product Type Botany 抗体

Research Area Botany

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted:Rice)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 52kDa

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from Oryza sativa Os05g0477200 protein:
221-320/472

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution (predicted:Rice)1M TBS(pH7.4) with 1% BSA, (predicted:Rice)3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com
