

Rabbit Anti-Shiga-like toxin IIe variant subunit A antibody

SL0880R

Product Name Shiga-like toxin IIe variant subunit A

Chinese Name 猪水肿病-志贺样毒素II型变异体 A 抗体

Alias Shiga-like toxin IIe variant subunit A; SLT-IIeA.

Research Area Bacteria and viruses

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted: Escherichia Coli,)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA

Applications (Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 35kDa

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from Shiga-like toxin IIe variant subunit A: 221-319/

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com
