

Rabbit Anti-Vibrio cholerae O1 serotype Ogawa antibody

SL0926R

Product Name	Vibrio cholerae O1 serotype Ogawa
Chinese Name	稻叶 O1 血清型霍乱弧菌抗体
Alias	Vibrio cholerae O1 serotype Ogawa. O1 群小川型霍乱弧菌; 霍乱弧菌;
Research Area	Bacteria and viruses
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Vibrio cholerae O1 serotype Ogawa) ELISA=1:5000-10000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration	1mg/ml
immunogen	whole cell protein of Vibrio cholerae O1 serotype Ogawa
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:Vibrio cholerae O1 serotype Ogawa)1M TBS(pH7.4) with 1% BSA, (predicted:Vibrio cholerae O1 serotype Ogawa)3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	