

## Rabbit Anti-Anthrax antibody

SL0935R

<b>Product Name</b>	Anthrax
<b>Chinese Name</b>	炭疽菌抗体
<b>Research Area</b>	Microbiology
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Bacillus anthraci) ELISA=1:5000-10000
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	whole cell protein of Bacillus anthraci
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	(predicted:Bacillus anthraci)1M TBS(pH7.4) with 1% BSA, (predicted:Bacillus anthraci)3% Proclin300 and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
<b>PubMed</b>	<a href="#">PubMed</a>
<b>Product Detail</b>	炭疽杆菌（Bacillus anthraci）属于需氧芽胞杆菌属，能引起羊、牛、马等动物及人类的炭疽病。