



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

Rabbit Anti-H7N9 Matrix Protein 2 antibody

SL10300R

Product Name H7N9 Matrix Protein 2

Chinese Name A 型禽流感病毒 H7N9 M1 蛋白抗体

Alias matrix protein 2 [Influenza A virus (A/northern shoverl/Mississippi/110S145/2011(H7N9))]; H7N9 M2.

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted:H7N9)

Applications WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 10kDa

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from H7N9 Matrix Protein 2: 51-97/97

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution (predicted:H7N9)1M TBS(pH7.4) with 1% BSA, (predicted:H7N9)3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail